

Features at a glance:

- for emergency, transport, clinic and home care
- for spot-check measurements
- compact, lightweight, simple and durable
- Masimo performance

sat 800

SpotCheck Pulse Oximeter

Made in Germany

Signal IQ™

rubber
protection

USB
service interface



easy-to-read
LC-display

simple and intuitive
operation



Package Contents

sat 800 O/N 36-1000

including:

- sat 800
- finger clip sensor with integrated cable
- instructions for use

accessories:

- external battery charger
- table stand

Specifications

dimensions:	128 x 85 x 46 mm
weight:	230 g incl. batteries
power supply:	2 AA batteries LR6
operating time:	20 hrs typ.
operating temperature:	+ 0° to + 50°C
storage temperature:	- 10° to + 70°C

SpO₂:	display:	1 - 100 %
	calibration:	70 - 100 %
	resolution:	1 %
pulse:	display:	25 - 240 bpm
	calibration:	25 - 240 bpm
	resolution:	1 bpm
quality:	Perfusion Index, Signal IQ™	
status:	sensor failure, low battery, system failure	

Masimo pulse oximetry works accurately where conventional pulse oximetry tends to fail. More than 110 peer-reviewed clinical studies document the significant superior performance of Masimo pulse oximetry, especially in challenging situations, such as motion and low perfusion.

version 2.3 - 12-2006 © Bitmos GmbH 2006

Bitmos GmbH
Himmelgeister Str. 37
D 40225 Düsseldorf, Germany
Phone: +49 211 60101030
Fax: +49 211 60101050
sales@bitmos.de
www.bitmos.de

Bitmos
MEDIZINTECHNIK