

Features at a glance:

- 30.000 operating hours warranty
- innovative maintenance concept
- simple operation
- incredibly silent

oxy 6000

Oxygen Concentrator

Made
in Germany



whisper silent
35 db(A)

modular
component
setup

microprocessor controlled

easily accessible
filters and fuses

USB interface

CE 0197

innovative multi-sensor

Scope of delivery

oxy 6000 O₂ concentrator

O/N 41-1000

including:

concentrator	41-1001
humidifier	41-1502
connector	41-1504
nasal cannula 2m	41-1506
nasal cannula 5m	41-1507
instructions for use	41-1201

more accessories on request

Specifications

dimensions:	535 x 203 x 520 mm
weight:	19.8 kg
power:	230 V / 50 Hz
power consumption:	350 W
fuses - external:	2 x TT 2.5 AH 250 V
fuses - internal:	T 1.0 AL 250 V
operating temperature:	+10° - +40°C
storage temperature:	-20° - +50°C
sound power:	35 dB(A)
MDD Classification:	IIa

flow: 0 to 2 lpm in 0.1 steps

2 to 4 lpm in 0.2 steps

4 to 6 lpm in 0.5 steps

oxygen: 1 to 4 lpm: 95 % - 3 %

4 to 5 lpm: 85 % +/- 3 %

5 to 6 lpm: 75 % +/- 3 %

alarms: temperature, power failure, flow,
oxygen, sensor, system

dust filter: in the rear

fine filter: behind maintenance lid

warranty: 30.000 hours or 5 years

Bitmos GmbH
Himmelgeister Str. 37
D 40225 Düsseldorf Germany

Phone: +49 211 60101030
Fax: +49 211 60101050

sales@bitmos.de
www.bitmos.de

Bitmos
MEDIZINTECHNIK